

HTL TIPS SPECIFICATION

Catalog number	Range of the tip volumes [μl]	Total tip length [mm]	Length of the handle part [mm]	External diameter of part cooperating with ejector [mm]	External diameter of part cooperating with box [mm]	Internal diameter of part cooperating with shaft (5,5mm away from the end of the tip) [mm]
301-03-052 - 301-03-152, 302-05-102	0,5-10	31,7	9,6	6,0	4,2	do 3
301-22-052 - 301-22-152, 302-26-152	0,5-10	45,7	12,3	5,9	4,2	3,3
302-03-102	0,5-20	50,4	16,0	7,3	4,7	5,0
302-08-102	1-100	51,0	16,0	7,3	4,8	4,9
301-06-052 - 301-79-052, 302-04-102	1-200	50,6	16,0	7,3	4,8	4,9
301-16-052, 301-16-102, 302-12-102	1-300	52,3	11,6	7,3	5,7	5,1
301-01-202 -301-01-352, 302-01-102	100-1000	78,0	20,2	9,3	7,9	7,3
301-71-052	1000-5000	123,5	21,0	14,9	13,0	12,4
301-75-052	1000-10000	143,8	21,0	17,7	15,9	15,1