



# One-Touch Vortexer

LX5611TV

[www.labdex.com](http://www.labdex.com)  
[info@labdex.com](mailto:info@labdex.com)

## One-Touch Vortexer LX5611TV

One-Touch Vortexer LX5611TV is a touch to mix operational unit suitable for mixing or centrifugal tubes. Equipped with a mixing speed variance of 0 to 4000 rpm and maximum sample capacity of 50 ml tube. Outer shell and touch concave are made up of chemically resistant plastic making this unit user friendly

### Features

- Maximum diameter of mixing tube - 30 mm
- Suitable for mixing and centrifugal tubes with mixing effect
- Chemically resistant outer shell and touch concave
- Speed variance is 0 ~ 4000 rpm
- Circular shaking orbit - 4.5 mm

### Applications

Used in molecular hybridization, blood banking and various samples requiring mixing

### Specifications

Model No	LX5611TV
Shaking orbit	4.5 mm ( circle )
Speed Range	0 ~ 4000 rpm
Maximum sample capacity	50 mL tube
Power input	1.2 W
Power output	0.8 W
Voltage	AC 100 ~ 230 V, 50 / 60 Hz
Power	12 W
Dimension	110 × 110 × 86 mm
Weight	1.2 kgs