



Centrifuge Rotors and Adapters



Centrifuges	01-800X		01-821X	01-830X				
Rotors	 							
Rotor Cat No.	01-8091	01-8092	01-8291	01-8393	01-8391	01-8394	01-8395	01-8392
Max. speed (rpm)	10000	10000	15000	15000	15000	15000	15000	15000
Max. RCF (×g)	5460	5460	15100	21380	21380	21380	21380	21380
Rotor capacity	2mL/1.5mL×8 0.2mL×8 0.5mL×8	0.2mL×16 PCR strips 0.2mL×2 PCR 8 strips	2mL/1.5mL×12 0.2mL×12 0.5mL×12	2mL/1.5mL×24 0.2mL×24 0.5mL×24	5mL×18	0.5mL×36 0.2mL×36	0.2mL/PCR8×4	5mL×12
Rotor material	High strength plastic	High strength plastic	High strength plastic	Aluminum alloy				
Bio-safe sealing				•	•	•	•	•
Rotor adapters	N/A				N/A		N/A	N/A
01-8093	•							
01-8094	•							
01-8293			•	•				
01-8294			•	•				
01-8396						•		

Centrifuge Rotors Adapters

Adapters	Cat No.	Description
	01-8093	2mL to 0.2mL rotor adapter, used with 01-8091 rotor, 8pcs/pk
	01-8094	2mL to 0.5mL rotor adapter, used with 01-8091 rotor, 8pcs/pk
		
	01-8396	0.5mL to 0.2mL rotor adapter, used with 01-8394 rotor, 36pcs/pk