



One-Touch Vortexer

LX5610TV

www.labdex.com
info@labdex.com

One-Touch Vortexer LX5610TV

One-Touch Vortexer LX5610TV is a touch to mix operational unit suitable for mixing or centrifugal tubes. Equipped with a mixing speed of 4000 rpm and maximum sample capacity of 50 ml tube. Outer shell and touch concave are made up of chemically resistant plastic making this unit user friendly

Features

- Maximum diameter of mixing tube - 30 mm
- Suitable for mixing and centrifugal tubes with mixing effect
- Chemically resistant outer shell and touch concave
- Circular shaking orbit - 4.5 mm

Applications

Used in molecular hybridization, blood banking and various samples requiring mixing

Specifications

Model No	LX5610TV
Shaking orbit	4.5 mm (circle)
Speed Range	0 ~ 4000 rpm
Maximum sample capacity	50 mL tube
Power input	1.2 W
Power output	0.8 W
Voltage	AC 100 ~ 230 V, 50 / 60 Hz
Power	12 W
Dimension	110 × 110 × 86 mm
Weight	1.2 kgs